

## XH-D Vortex Mixer



### Features

Adapted to re-suspension particles, extract drugs and vortex cell suspensions. Mixing tissue samples. And also can be used to ELISA reaction and RIA Tests.

Multiple plates available, can be used for various kinds of test tubes and small containers.

High performance & low noise AC motor to make rapid and thorough mixing.

Continuous and Inching model available, Easy to operate, Knob design makes it easy carry out stepless speed regulating.

Reasonable design, compact structure and stable performance.

### Main Technical Data

Operation method	Circling
Circling diameter	4.4mm
Rating speed	1000-2800rpm
Running mode	Constant, Inching
Working plate size	—
Input voltage(frequency)	AC 220V±10%(50/60Hz)
Permissible relative moisture	80%
Output power	30W
Protective class	IP21
Permissible ambient temperature:	5-50℃
Dimensions	180×120×180mm
Net. Weight	4.5Kg